


Product Datasheet

Rat Galectin-1 ELISA Kit (orb507516)

Description	Rat Galectin-1 ELISA Kit
Reactivity	Rat
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A2-E135
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is 780-3000 pg/ml
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	P11762
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Rat Galectin-1 PicoKine ELISA Kit standa...