


Product Datasheet

Human Serpin A1/alpha 1-Antitrypsin ELISA Kit (orb507524)

Description	Human Serpin A1/alpha 1-Antitrypsin ELISA Kit
Reactivity	Human
Range	1.56 ng/ml - 100 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E25-K418
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:10000, concentration in serum and plasma is 22.2-270 ug/ml
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and urine.
Uniprot ID	P01009
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human Serpin A1/alpha 1-Antitrypsin Pico...