


## Product Datasheet

### Human TRKC ELISA Kit (orb507525)

<b>Description</b>	Human TRKC ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	156 pg/ml - 10,000 pg/ml
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Expression system for standard: NS0; Immunogen sequence: C32-A428
<b>Storage</b>	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
<b>Note</b>	For research use only
<b>Application notes</b>	Dilution ratio of 1:1, concentration in serum and plasma is 1-1.6 ng/ml
<b>Sample Types</b>	cell culture supernatants, serum and plasma (heparin, EDTA).
<b>Uniprot ID</b>	<b>Q16288</b>
<b>Sensitivity</b>	10 pg/ml
<b>Expiration Date</b>	Please enquire.



Human TRKC  
PicoKine  
ELISA Kit  
Standard C...