


Product Datasheet

Human Transthyretin / TTR / PreAlbumin ELISA Kit (orb507527)

| | |
|----------------------------|---|
| Description | Human Transthyretin / TTR / PreAlbumin ELISA Kit |
| Reactivity | Human |
| Range | 312 pg/ml - 20,000 pg/ml |
| Tested Applications | ELISA |
| Immunogen | Expression system for standard: E.coli; Immunogen sequence: G21-E147 |
| Storage | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| Note | For research use only |
| Application notes | Dilution ratio of 1:1000000, concentration in serum and plasma is around 240 ug/ml |
| Sample Types | cell culture supernatants, serum and plasma (heparin, EDTA). |
| Uniprot ID | P02766 |
| Sensitivity | 10 pg/ml |
| Expiration Date | Please enquire. |



Human Transthyretin/TTR PicoKine ELISA K...