


Product Datasheet

Rat PDGF-BB ELISA Kit (orb507539)

Description	Rat PDGF-BB ELISA Kit
Reactivity	Rat
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: S82-T190
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	For serum, dilution ratio is 1:1 and concentration is 1643-4291 pg/ml. For plasma, dilution ratio is 1:1 and concentration is 174-1861 pg/ml
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q05028
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Rat PDGF-BB
PicoKine
ELISA Kit
standard ...