


Product Datasheet

Human FGFBP1 ELISA Kit (orb507544)

Description	Human FGFBP1 ELISA Kit
Reactivity	Human
Range	23.4 pg/ml - 1500 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: K24-C234
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:10, concentration in serum and plasma is around 1.7 ng/ml
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q14512
Sensitivity	5 pg/ml
Expiration Date	Please enquire.



Human
FGFBP1
PicoKine
ELISA Kit
standard...