

Product Datasheet

DOK1 (phospho-Tyr398) antibody (orb5085)

Description	DOK1 (phospho-Tyr398) antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human DOK1 around the phosphorylation site of Tyr398 (EG(-Y)EL)
Target	phospho-DOK1 (Tyr398)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	53 kDa
Uniprot ID	Q99704
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.