

Product Datasheet

RPA1 antibody (Biotin) (orb51009)

| | |
|----------------------------|--|
| Description | RPA1 antibody (Biotin) |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Biotin |
| Tested Applications | ELISA |
| Immunogen | Recombinant Human Replication protein A 70 kDa DNA-binding subunit protein (2-616AA) |
| Target | RPA1 |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P27694 |
| Dilution Range | ELISA |
| Expiration Date | 12 months from date of receipt. |