


Product Datasheet

Chicken Ovotransferrin protein (orb510858)

Description	Recombinant chicken Ovotransferrin protein
Tested Applications	ELISA, SDS-PAGE, WB
Target	Chicken Ovotransferrin
Preservatives	PBS, pH 7.4
Form/Appearance	Liquid
Storage	Aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	Tag Info: Native protein Sequence Info: Full length protein
Purity	> 95%(SDS-PAGE)
Source	Purified from chicken egg white
Hazard Information	For research use only
Expiration Date	6 months from date of receipt.



(Tris-Glycine gel)
Discontinuous SDS-PAG...