



Product Datasheet

HIF-1alpha antibody (orb511147)


Description	HIF-1alpha antibody
Species/Host	Rabbit
Reactivity	Guinea pig, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human HIF-1alpha . Please contact us for the exact immunogen sequence.
Target	HIF-1alpha
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
Uniprot ID	Q16665
Dilution Range	IHC-P: 1:100-500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



IHC-P image of rat lymph node tissue usi...



IHC-P image of guinea pig bladder tissue...