



Product Datasheet

TYMS antibody (orb51122)


Description	TYMS antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IP, WB
Immunogen	Recombinant Human Thymidylate synthase protein (2-54AA)
Target	TYMS
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04818
Dilution Range	WB:1:500-5000, IF:1:50-200, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2(lane 1), ...



Immunofluorescence analysis of HeLa cell...



Immunoprecipitation analysis of HeLa who...