


Product Datasheet

Human Flt-3 Ligand Protein (orb511237)

Description	Human Flt-3 Ligand, His Tag is expressed from Baculovirus-Insect cells. It contains AA Thr 27 - Pro ...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Flt-3 Ligand
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 19.9 kDa. The protein migrates as 22-28 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	90%
Source	Human Flt-3 Ligand, His Tag (orb511237) is expressed from Baculovirus-Insect cells. It contains AA Thr 27 - Pro 185 (Accession # P49771-1).
MW	19.9 kDa
Uniprot ID	P49771
NCBI	NP_001450.2
Expiration Date	6 months from date of receipt.



Immobilized
Human Flt-3,
Fc Tag (Cat.
No...