

## Product Datasheet

### Mouse Siglec-2 / CD22 Protein (orb511261)

<b>Description</b>	Mouse Siglec-2, His Tag is expressed from human 293 cells (HEK293). It contains AA Ser 22 - Arg 702 ...
<b>Reactivity</b>	Mouse
<b>Endotoxins</b>	1.0 EU per µg
<b>Conjugation</b>	Unconjugated
<b>Target</b>	Siglec-2
<b>Preservatives</b>	PBS, pH7.4
<b>Form/Appearance</b>	Powder
<b>Storage</b>	-20°C
<b>Tag</b>	C-10×His
<b>Note</b>	For research use only
<b>Application notes</b>	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 78.3 kDa. The protein migrates as 100 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
<b>Protein Sequence</b>	P35329.1
<b>Purity</b>	95%
<b>Source</b>	Mouse Siglec-2, His Tag (orb511261) is expressed from human 293 cells (HEK293). It contains AA Ser 22 - Arg 702 (Accession # P35329-1).
<b>MW</b>	78.3 kDa
<b>Uniprot ID</b>	<a href="#">P35329</a>
<b>NCBI</b>	<a href="#">P35329.1</a>
<b>Expiration Date</b>	6 months from date of receipt.