

Product Datasheet

Mouse PLGF / PGF Protein (orb511267)

Description	Mouse PLGF, His Tag is expressed from human 293 cells (HEK293). It contains AA Val 19 - Pro 158 (Acc...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	PLGF
Preservatives	100 mM HAC, pH2.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 17.8 kDa. The protein migrates as 27-35 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_032853.1
Purity	95%
Source	Mouse PLGF, His Tag (orb511267) is expressed from human 293 cells (HEK293). It contains AA Val 19 - Pro 158 (Accession # P49764-1).
MW	17.8 kDa
Uniprot ID	P49764
NCBI	NP_032853.1
Expiration Date	6 months from date of receipt.