
Product Datasheet

Human TSLP Protein (orb511287)

Description

Human TSLP Protein, His Tag, premium grade

It contains A...

Reactivity

Human

Endotoxins

0.1 EU per µg

Conjugation

Unconjugated

Target

TSLP

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 16.8 kDa. The protein migrates as 20-27 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under non-reducing (NR) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

NP_149024.1

Purity

90%

Source

Human TSLP, His Tag (orb511287) is expressed from human 293 cells (HEK293). It contains AA Tyr 29 - Gln 159 (Accession # Q969D9-1).

MW

16.8 kDa

Uniprot ID[Q969D9](#)**NCBI**[NP_149024.1](#)**Expiration Date**

6 months from date of receipt.