

Product Datasheet

Rhesus macaque PLGF / PGF Protein (orb511301)

Description	Rhesus macaque PLGF, His Tag is expressed from human 293 cells (HEK293). It contains AA Leu 19 - Arg...
Reactivity	Monkey
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	PLGF
Preservatives	100 mM Acetic Acid, pH2.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 19.3 kDa. The protein migrates as 28-34 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	XP_001096010.1
Purity	90%
Source	Rhesus macaque PLGF, His Tag (orb511301) is expressed from human 293 cells (HEK293). It contains AA Leu 19 - Arg 170 (Accession # F7HB10-1).
MW	19.3 kDa
Uniprot ID	F7HB10
NCBI	XP_001096010.1
Expiration Date	6 months from date of receipt.