

---

## Product Datasheet

### Cynomolgus LILRB4 / CD85k / ILT3 Protein (orb511303)

**Description**

Cynomolgus LILRB4, Mouse IgG2a Fc

(HEK293). It contains AA Gln...

<b>Reactivity</b>	Monkey
<b>Endotoxins</b>	1.0 EU per µg
<b>Conjugation</b>	Unconjugated
<b>Target</b>	LILRB4
<b>Preservatives</b>	50 mM Tris, 100 mM Glycine, 150 mM NaCl, pH7.5
<b>Form/Appearance</b>	Powder
<b>Storage</b>	-20°C
<b>Tag</b>	C-mFc(mIgG2a)
<b>Note</b>	For research use only
<b>Application notes</b>	This protein carries a mouse IgG2a Fc tag at the C-terminus. The protein has a calculated MW of 53.2 kDa. The protein migrates as 65 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
<b>Protein Sequence</b>	XP_015297198.1
<b>Purity</b>	95%
<b>Source</b>	Cynomolgus LILRB4, Mouse IgG2a Fc Tag (orb511303) is expressed from human 293 cells (HEK293). It contains AA Gln 22 - Glu 259 (Accession # XP_015297198.1).
<b>MW</b>	53.2 kDa
<b>Uniprot ID</b>	XP_015297198.1
<b>NCBI</b>	<a href="#">XP_015297198.1</a>
<b>Expiration Date</b>	6 months from date of receipt.