

Product Datasheet

Human OX40 / TNFRSF4 / CD134 Protein (orb511311)

Description

Human OX40 Protein, Mouse IgG2a Fc

(HEK293). It contains AA Le...

Reactivity

Human

Endotoxins

0.1 EU per µg

Conjugation

Unconjugated

Target

OX40

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-mFc(mlgG2a)

Note

For research use only

Application notes

This protein carries a mouse IgG2a Fc tag at the C-terminus. The protein has a calculated MW of 47.1 kDa. The protein migrates as 58-66 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

NP 003318.1

Purity

95%

Source

Human OX40 Protein, Mouse IgG2a Fc Tag (orb511311) is expressed from human 293 cells (HEK293). It contains AA Leu 29 - Ala 216 (Accession # P43489-1).

MW

47.1 kDa

Uniprot ID
P43489
NCBI
003318.1
Expiration Date

6 months from date of receipt.