

---

## Product Datasheet

### Human TSLP R Protein (orb511319)

**Description**

Human TSLP R, Fc Tag is expressed from

25 - Lys 231 (Ac...

**Reactivity**

Human

**Endotoxins**

1.0 EU per µg

**Conjugation**

Unconjugated

**Target**

TSLP R

**Preservatives**

Tris with Glycine, Arginine and NaCl, pH7.5

**Form/Appearance**

Powder

**Storage**

-20°C

**Tag**

C-Fc

**Note**

For research use only

**Application notes**

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 50.4 kDa. The protein migrates as 60-65 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Protein Sequence**

NP\_071431.2

**Purity**

95%

**Source**

Human TSLP R, Fc Tag (orb511319) is expressed from human 293 cells (HEK293). It contains AA Gly 25 - Lys 231 (Accession # Q9HC73-1).

**MW**

50.4 kDa

**Uniprot ID**

[Q9HC73](#)

**NCBI**

[NP\\_071431.2](#)

**Expiration Date**

6 months from date of receipt.