

Product Datasheet

PIK3CG antibody (Biotin) (orb51213)

| | |
|----------------------------|--|
| Description | PIK3CG antibody (Biotin) |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Biotin |
| Tested Applications | ELISA |
| Immunogen | Recombinant Human Phosphatidylinositol 4,5-bisphosphate 3-kinase catalytic subunit gamma isoform protein (438-588AA) |
| Target | PIK3CG |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P48736 |
| Dilution Range | ELISA |
| Expiration Date | 12 months from date of receipt. |