

Product Datasheet

OLE9 antibody (HRP) (orb51263)

Description	OLE9 antibody (HRP)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Olea europaea Glucan endo-1,3-beta-D-glucosidase protein (360-460AA)
Target	OLE9
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q94G86
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.