


Product Datasheet

EGFR (phospho-Tyr1197) antibody (orb5128)

Description	EGFR (phospho-Tyr1197) antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human EGFR around the phosphorylation site of Tyr1197 (AE(-p-Y)LR)
Target	phospho-EGFR (Tyr1197)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	170 kDa
Uniprot ID	P00533
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat live...