

Product Datasheet

GP antibody (HRP) (orb51291)

Description	GP antibody (HRP)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Zaire ebolavirus Pre-small/secreted glycoprotein protein (33-364AA)
Target	GP
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P87670
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.