

Product Datasheet

Defensin-like protein 1 antibody (Biotin) (orb51309)

Description	Defensin-like protein 1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant <i>Dahlia merckii</i> Defensin-like protein 1 protein (1-50AA)
Target	Defensin-like protein 1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0C8Y4
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.