

Product Datasheet

Sporamin A antibody (FITC) (orb51324)

Description	Sporamin A antibody (FITC)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	FITC
Immunogen	Recombinant Ipomoea batatas Sporamin A protein (24-219AA)
Target	Sporamin A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P14715
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.