

Product Datasheet

TTC30A antibody (orb514155)

Description	Rabbit polyclonal antibody to TTC30A
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant fragment corresponding to 272-466aa of human IFT70A protein
Target	IFT70A
Form/Appearance	Lyophilized powder with 1% rAlbumin and 0.02% NaN ₃ .
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Reconstitute in 100µL of sterile H ₂ O with 50% glycerol.
Isotype	IgG
Clonality	Polyclonal
MW	76kDa
Dilution Range	WB: 1:200-1:1,000
Expiration Date	12 months from date of receipt.