

Product Datasheet

NDI1 antibody (orb515166)

Description	Rabbit polyclonal antibody to NDI1
Species/Host	Rabbit
Reactivity	Yeast
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	A peptide corresponding to the C terminal of <i>Saccharomyces cerevisiae</i> NDI1 protein
Target	NDI1
Form/Appearance	Lyophilized powder with 0.02% NaN ₃ .
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Reconstitute in 100µL of sterile H ₂ O with 50% glycerol.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	ELISA: 1:20,000-1:80,000
Expiration Date	12 months from date of receipt.