

Product Datasheet

TIE1 antibody (orb515176)

Description	Mouse monoclonal antibody to TIE1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	A peptide corresponding to middle region of human TIE1 protein
Target	TIE1
Form/Appearance	Lyophilized powder with 0.02% NaN ₃ .
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Reconstitute in 100µL of sterile H ₂ O with 50% glycerol.
Isotype	IgG2a
Clonality	Monoclonal
MW	60kDa
Dilution Range	WB: 1:200-1:1,000
Expiration Date	12 months from date of receipt.