



## Product Datasheet

### ERN1 antibody (orb515588)

<b>Description</b>	Mouse monoclonal antibody to ERN1
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IH, WB
<b>Immunogen</b>	Purified recombinant fragment of human IRE1 expressed in E. Coli.
<b>Target</b>	ERN1
<b>Preservatives</b>	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
<b>Form/Appearance</b>	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>Source</b>	Mouse
<b>Uniprot ID</b>	<a href="#">O75460</a>
<b>Entrez</b>	<a href="#">2081</a>
<b>Dilution Range</b>	WB: 1:500-1000, IHC-P: 1:100-500
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of IRE1 expression...



Immunohistochemical analysis of IRE1 sta...