

Product Datasheet

FCGR3B antibody (Biotin) (orb515677)

Description	FCGR3B antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Low affinity immunoglobulin gamma Fc region receptor III-B protein (18-125AA)
Target	FCGR3B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O75015
Expiration Date	12 months from date of receipt.