

Product Datasheet

MAP4K2 antibody (FITC) (orb515724)

Description	MAP4K2 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Mitogen-activated protein kinase kinase kinase 2 protein (386-469AA)
Target	MAP4K2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q12851
Expiration Date	12 months from date of receipt.