

Product Datasheet

DLGAP2 antibody (FITC) (orb516207)

Description	DLGAP2 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Disks large-associated protein 2 protein (31-139AA)
Target	DLGAP2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9P1A6
Expiration Date	12 months from date of receipt.