

Product Datasheet

IBA57 antibody (Biotin) (orb516704)

Description	IBA57 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Putative transferase CAF17, mitochondrial protein (286-356AA)
Target	IBA57
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5T440
Expiration Date	12 months from date of receipt.