

Product Datasheet

MEX3A antibody (Biotin) (orb517254)

Description	MEX3A antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human RNA-binding protein MEX3A protein (416-466AA)
Target	MEX3A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A1L020
Expiration Date	12 months from date of receipt.