

Product Datasheet

rbm47 antibody (HRP) (orb517360)

Description	rbm47 antibody (HRP)
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Zebrafish RNA-binding motif protein 47 protein (390-524AA)
Target	rbm47
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A2BIC5
Expiration Date	12 months from date of receipt.