

Product Datasheet

rhoac antibody (Biotin) (orb517719)

Description	rhoac antibody (Biotin)
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Zebrafish Rho-related GTP-binding protein RhoA-C protein (1-190AA)
Target	rhoac
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7T2E8
Expiration Date	12 months from date of receipt.