


Product Datasheet

mycn antibody (orb517808)

Description	mycn antibody
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Zebrafish N-myc protein (1-93AA)
Target	mycn
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9PSJ0
Expiration Date	12 months from date of receipt.



Western blot
analysis of
Recombinant
pr...