

Product Datasheet

ARNTL antibody (FITC) (orb517858)

| | |
|------------------------|---|
| Description | ARNTL antibody (FITC) |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | FITC |
| Immunogen | Recombinant Human Aryl hydrocarbon receptor nuclear translocator-like protein 1 protein (366-585AA) |
| Target | ARNTL |
| Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | O00327 |
| Expiration Date | 12 months from date of receipt. |