


Product Datasheet

GNAT2 antibody (orb517888)

Description	GNAT2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Guanine nucleotide-binding protein G(t) subunit alpha-2 protein (2-139AA)
Target	GNAT2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P19087
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of HeLa cell...