


Product Datasheet

HOXC6 antibody (orb517896)

Description	HOXC6 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Homeobox protein Hox-C6 protein (1-120AA)
Target	HOXC6
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P09630
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of PC3 cells...