

Product Datasheet

ORF7 antibody (HRP) (orb518105)

Description	ORF7 antibody (HRP)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Porcine reproductive and respiratory syndrome virus Nucleocapsid protein (1-123AA)
Target	ORF7
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	V5N6H2
Expiration Date	12 months from date of receipt.