

Product Datasheet

BLLF3 antibody (Biotin) (orb518159)

Description	BLLF3 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Epstein-Barr virus Deoxyuridine 5\'-triphosphate nucleotidohydrolase protein (1-278AA)
Target	BLLF3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	K9US42
Expiration Date	12 months from date of receipt.