


Product Datasheet

NS antibody (orb518212)

Description	NS antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Influenza A virus Non-structural protein 1 protein (1-237AA)
Target	NS
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q0HD54
Expiration Date	12 months from date of receipt.



Western blot analysis of Positive WB de...