

Product Datasheet

P antibody (Biotin) (orb518223)

Description	P antibody (Biotin)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Vesicular stomatitis Indiana virus Protein C\ protein (1-67AA)
Target	P
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q86132
Expiration Date	12 months from date of receipt.