

Product Datasheet

Parasite GST protein (orb5186)

| | |
|------------------------|--|
| Description | Recombinant Schistosoma Japonicum Glutathione S-transferase Class-mu 26 kDa Isozyme is produced by o... |
| Endotoxins | Less than 0.1 ng/ug (1 EU/ug) |
| Form/Appearance | Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4. |
| Storage | Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months. |
| Tag | tag free |
| Note | For research use only |
| Uniprot ID | P08515 |
| Expiration Date | 6 months from date of receipt. |