

Product Datasheet

Ferritin-2, chloroplastic antibody (FITC) (orb51868)

Description	Ferritin-2, chloroplastic antibody (FITC)
Species/Host	Rabbit
Reactivity	Plant
Conjugation	FITC
Immunogen	Recombinant Glycine max Ferritin-2, chloroplastic protein (52-257AA)
Target	Ferritin-2, chloroplastic
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q94IC4
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.