


Product Datasheet

PFN1 antibody (orb519391)

| | |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Description | PFN1 antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, WB |
| Immunogen | Fusion protein of Human PFN1 |
| Target | PFN1 |
| Preservatives | -20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P07737 |
| Dilution Range | ELISA:1:1000-1:2000, WB:1:200-1:1000 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of Lane 1-3: HeLa...