

Product Datasheet

PTPN11 Antibody (orb519476)

Catalog Number	orb519476
Category	Antibodies
Description	This PTPN11 Antibody is an unconjugated polyclonal product. It targets PTPN11 using a Fusion protein of Human PTPN11 as the immunogen. This antibody is suitable for ELISA, IHC, WB. Purification: Antigen affinity purification.
Target	PTPN11
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Human, Mouse, Rat
Form/Appearance	Liquid
Buffer/Preservatives	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Purification	Antigen affinity purification
Immunogen	Fusion protein of Human PTPN11
UniProt ID	Q06124
Tested applications	ELISA, IHC, WB
Dilution range	ELISA:1:1000-1:2000, WB:1:200-1:1000, IHC-P:1:25-1:100

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States


Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note


For research use only



Immunohistochemistry of paraffin-embedded Human gastric cancer tissue using PTPN11 antibody



Western blot aanalysis of Lane 1-2: NIH/3T3 cells, hela cells using PTPN11 antibody



Immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using PTPN11 antibody

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558