


Product Datasheet

FAK (phospho-Tyr407) antibody (orb5213)

Description	FAK (phospho-Tyr407) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human FAK around the phosphorylation site of Tyr407 (DT(p-Y)TM)
Target	Phospho-FAK (Tyr407)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	116 kDa
Uniprot ID	Q05397
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, Flow-Cyt=1µg /Test, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Flow cytometry of Hep G2(blue line).Prim...