



Product Datasheet

TRIM39 antibody (orb52142)

Description	TRIM39 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human E3 ubiquitin-protein ligase TRIM39 protein (310-416AA)
Target	TRIM39
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9HCM9
Dilution Range	WB:1:500-5000, IHC-P:1:500-1000, IF:1:50-500
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Immunofluorescent analysis of HeLa cells...